

9.3 Error code / Error chart / Remedy

| | Error / Cause | Display |
|------------|--|--|
| 100 | <p>On door US, Reed sensor S 1 or auxiliary relays K1 is not closed or the potentiometer value is not within the proper range or has changed during the current program.</p> <p>Query input S110/111 and R1/2 only during current program.</p> <p>Cause:</p> <p>Door switch S1 is closed S110</p> <p>Signal cable linear drive interrupted</p> <p>Linear drive</p> <p>Relay K1</p> <p>Actual value pot meters not in target range</p> <p>Chack all plug connections</p> <p>Remedy:</p> <p>defect, cable break, adjust</p> <p>interrupt of cable of connector</p> <p>replace</p> <p>replace</p> <p>readjust doors</p> | <p>Door is open</p> <p>Error code 100</p> |
| 101 | <p>(Only for WD 250/) The opening sequence of US door takes more than 20 seconds or the linear drive is defective.</p> <p>Control system will interrupt the opening sequence and display error if either the opening process for door with lifting motor US (M 35) takes more than 16 seconds until the lower Reed S 112) or the potentiometer value remains within the 'closed' range.</p> <p>Query input S 112, R1 only during opening procedure.</p> <p>Cause:</p> <p>Manual linear drive on door panel print (buttons)</p> <p>Open</p> <p>Activate door drive motor M 35 with control K35</p> <p>Reed sensor S 112 does not open</p> <p>Check all plug connections</p> <p>Safety relay ESR 1</p> <p>Relay ESR 1 LED green</p> <p>LED yellow = resistance incorrect</p> <p>LED red = cable break</p> <p>Remedy:</p> <p>re-calibrate door, exchange</p> <p>Linear motor</p> <p>replace door motor or capacitor</p> <p>replace S 112</p> <p>re-crimp S 112 to print</p> <p>no action required</p> <p>replace safety strip</p> <p>replace safety strip or find</p> <p>Cable break</p> | <p>Door does not open</p> <p>Error code 101</p> |
| 102 | <p>The closing process of the US door takes more than 20 seconds or linear drive is defective.</p> <p>The control system will interrupt closing sequence and displays error if either the closing procedure of door with lifting motor US (M 35) takes in excess of 16 seconds until the upper Reed sensor (S 1/K1) closes or until the potentiometer value R1 comes within 'closed' range.</p> <p>Query of input S 110, R1 only during closing process.</p> | <p>Door does not close</p> <p>Error code 102</p> |
| 103 | <p>The safety strip of door US has triggered.</p> <p>The control system will interrupt closing sequence and displays error if the closing procedure of door with lifting motor US (M 35) the safety strip (S 209) opens the relay ESR1 and reopens the</p> | <p>Object in door</p> <p>Error code 103</p> |

| | | |
|------------|--|---|
| | door (M 35). Query input S109 only during closing procedure. | |
| 105 | On door US , the Reed sensor or the auxiliary relay S 2 is not closed or the UCvalue is not within the 'closed' range or has opened or has gone outside of range during the current program. Program will shut down under error display. Query input S111, R2 only during current program. In case of single door model (Pos. 5.11), there will be no query. | Door is open Error code 105 |
| 106 | The opening process of door CS takes more than 20 seconds. If opening process of door with lifting motor US (M 37) takes more than 16 seconds until the lower Reed sensor (S 113) opens, or if potentiometer R2 value does not go out of 'closed' range, the control system will interrupt the opening process under display of error. Query input S113, R2 only during opening process. In case of single door version (Pos. 5.11), there is no query. | Door does not open Error code 106 |
| 107 | The closing process of door CS takes more than 20 seconds. If closing process of door with lifting motor CS (M 37) takes more than 16 seconds until the upper Reed sensor (S2/K2) opens, or if potentiometer R2 does not come into 'closed' range, the control system interrupts the closing process and displays error. Query input S111 only during closing process. | Door does not close Error code 107 |
| 108 | The safety strip of door RS has triggered. If during the closing process of door CS (M 37), the safety strip (S 210) opens the relay ESR2, the control system will interrupt the closing process and display error and reopens the door (M 37). Query input S110 only during closing process. Cause: Relay ESR 1 green LED yellow = resistance incorrect LED red = cable break Safety relay ESR 1 is defective Remedy: no action required replace safety strip replace safety strip or find cable break replace ESR 1 | Object in door Error code 108 |
| 109 | A communications error between processor and door press prints has occurred or the door press print has a defect. Cause: Door pressure print defective Fuse blown Connector plug has loose contact Coding switch sitting in wrong position Remedy Replace door pressure print Check connectors, fuse and coding switch and correct -> 8 | No communication Error code 109 |
| 110 | The minimal pressure in the wash system has been undercut. If the pressure switch 150 m bar (S 211) designed to monitor pump pressure, remains open for 90 seconds past start of fill phase, the control system will interrupt the current program and display error. Query input S211 only active with media M1-4, M7-9 for 90 seconds after program time and under 85 °C wash liquid temperature. Over 85°C, query is inactive. | |

| | | |
|------------|---|--|
| | <p>Cause:</p> <p>Rotation direction of wash pump</p> <p>Foam:</p> <p>Incorrect conn. of wash additive</p> <p>Pump pressure switch S 211 does not close</p> <p>Remedy:</p> <p>Reverse rotation direction by exch. Phase 1 and 2</p> <p>flush chamber with warm water</p> <p>correct container connection</p> <p>replace</p> | |
| 112 | <p>Floor pan has collected water.</p> <p>If input SB (optional float switch) is longer then 2 seconds closed, the control system will interrupt the current program and display error.</p> <p>Query input SB also in "Program ready" mode.</p> <p>Cause:</p> <p>Coarse or fine mesh sieve clogged</p> <p>Circuit over switch S B is open</p> <p>Leakage in system</p> <p>Remedy:</p> <p>clean sieves</p> <p>locate interruption</p> <p>fix leakage, dry basin</p> | <p>Machine leaks</p> <p>Error code 112</p> |
| 113 | <p>If input S 109 (opener) pressure switch opens - at 35 mm water column – within 2 seconds, display shows malfunction "Sieve clogged", current program is interrupted</p> <p>Query input differential pressure switch (normally closed) S 109 during entire program cycle with water M1-M4 and M9 only active during pre-rinse cycle F1.</p> <p>If an error occurs during pre-rinse, the program automatically repeats the pre-rinse cycle without program interruption. If the error occurs a second time, the current program is interrupted by means of 113. See function F1.</p> <p>Cause:</p> <p>Coarse or fine mesh sieve clogged</p> <p>Circuit over switch S 109 is open</p> <p>foam is being formed</p> <p>In case of significant loading with blood, the machine detects foam formation and automatically directs a repeat of the pre-rinse cycle. High quality pre-rinsing is an important prerequisite for efficient cleaning.</p> <p>Note:</p> <p>Foam is being formed that has an adverse effect on cleaning efficiency. Blood remnants left after the pre-rinse cycle enhance so-called soaping during the subsequent alkaline cleaning cycle. That results in reduced cleaning efficiency which increases the prion risk.</p> <p>Remedy:</p> <p>clean sieves</p> <p>locate interruption</p> <p>Change programparameter</p> | <p>Screen clogged</p> <p>Error code 113</p> |
| 121 | <p>The level sensor is defective or not connected.</p> <p>If level sensor (P 1) is defective or not connected for more than 2 seconds (the value must be between 0.2 and 4.5 V), the control system interrupts the current program and displays error.</p> <p>Query limit values input P1 during program cycle.</p> | <p>No water</p> <p>Error code 121</p> |
| 130 | <p>The temperature sensor on bottom of wash chamber is defective or not connected.</p> <p>If the temperature sensor NTC 1 For the wash chamber is defective or not connected (the value must be between 0 and 150 °C), the control system interrupts the current program and displays error.</p> | <p>NTC 1 defective</p> <p>Error code 130</p> |

| | | |
|------------|---|---|
| | <p>Query of limit value input NTC 1 during program cycle.</p> <p>Cause:</p> <p>If analysis T 1 = sensor defective</p> <p>Resistance readings:</p> <p>20°C = 12.490 K Ohm, 25°C = 10.000 K Ohm, 30°C = 8.060 K Ohm, 35°C = 6.356 K Ohm</p> <p>40°C = 5.331 k Ohm, 45°C = 4.373 K Ohm, 50°C = 3,606 K Ohm, 55°C = 2.989 K Ohm</p> <p>60°C = 2.949 K Ohm, 65.°C = 2.085 K Ohm, 70°C = 1.753 K Ohm, 75°C = 1.481 K Ohm</p> | |
| 131 | <p>The temperature sensor in the dryer is defective or not connected.</p> <p>See malfunction 130</p> | <p>NTC 2 defective</p> <p>Error code 131</p> |
| 132 | <p>The temperature sensor NTC 3 is defective or not connected.</p> <p>See malfunction 130</p> | <p>NTC 3 defective</p> <p>Error code 132</p> |
| 133 | <p>The temperature sensor in the DI boiler is defective or not connected.</p> <p>See malfunction 130</p> | <p>NTC 4 defective</p> <p>Error code 133</p> |
| 134 | <p>The Temperature sensor on bottom of wash chamber is defective or not connected.</p> <p>See malfunction 130</p> | <p>NTC 5 defective</p> <p>Error code 134</p> |
| 135 | <p>The temperature sensor in the DI boiler is defective or not connected.</p> <p>See malfunction 130</p> | <p>NTC 6 defective</p> <p>Error code 135</p> |
| 136 | <p>The temperature sensor on bottom of wash chamber is defective or not connected.</p> <p>If temperature sensor PT100 1 is defective or not connected (value must be between 0 and 150 °C), the control system will interrupt the current program and display error.</p> <p>Monitoring limit value input PT100 during program cycle.</p> | <p>PT100 1 defective</p> <p>Error code 136</p> |
| 137 | <p>Temperature sensor in dryer is defective or not connected.</p> | <p>PT100 2 defective</p> <p>Error code 137</p> |
| 138 | <p>Temperature sensor PT100 is defective or not connected.</p> | <p>PT100 3 defective</p> <p>Error code 138</p> |
| 139 | <p>Temperature sensor PT100 3 is defective or not connected.</p> | <p>PT100 4 defective</p> <p>Error code 139</p> |
| 140 | <p>Float switch SD in DI boiler remains open even 30 minutes following start of filling process. No DI water intake to boiler.</p> <p>Query input S SD only during intake process.</p> <p>Cause:</p> <p>DI water valve closed</p> <p>DI valve defective</p> <p>Float switch SD is defective</p> | <p>No DI water</p> <p>Error code 140</p> <p>Remedy:</p> <p>open DI water valve</p> <p>clean or replace</p> <p>replace SD (de-install heater)</p> |
| 140 | <p>Adjustment of limit value for query input S SD during intake process.</p> | <p>DI boiler SD</p> |

| | | |
|------------|--|---|
| | 10 – 45 min., default = 15 min. | 15 Minutes |
| 141 | <p>Working level with CW during intake process has not been attained. No cold water intake into wash chamber or door contact of door (S1 or S2) is open.</p> <p>If during filling of wash chamber with CW (valve Y 11) the operating level 2 of level sensor P1 is not attained within a predetermined period of time tcw (5 min.), the control system will interrupt the current program and display error.</p> <p>Monitoring limit value: only during intake process until operating level 2 is attained.</p> <p>Cause:</p> <p>Activate with Y 11</p> <p>CW valve closed</p> <p>CW valve defective</p> <p>Level monitoring system P1 leaks</p> <p>S 110 / S 111 door switch remains open</p> <p>Remedy:</p> <p>open CW valve</p> <p>clean replace valve</p> <p>check hoses</p> <p>check door switch and position</p> | No CW water Error code 141 |
| 141 | <p>Adjustment of limit value for cold water intake.</p> <p>1 - 15 min., default = 5 min.</p> | CW water P1 5 Minutes |
| 142 | <p>During WW intake process, the operating level is not attained. No warm water intake into wash chamber or door contact of door (S1 or S2) is open.</p> <p>If during WW intake into wash chamber (valve Y 12), working level 2 of level sensor P1 is not attained within a predetermined period of time tww (5 min.), the control system will interrupt the current program and display error.</p> <p>Monitoring limit value: As under 141</p> <p>See 141</p> | No WW water Error code 142 |
| 142 | <p>Adjustment of limit value for warm water intake tww.</p> <p>1 - 15 min., Default = 5 min.</p> | WW water P1 5 Minutes |
| 143 | <p>During intake process with DI water, working level is not attained.</p> <p>If during intake process with DI water (valve Y 13), operating level 2 is not attained within a predetermined period of time tvek (7 min.) the control system will interrupt the current process and display error.</p> <p>Monitoring limit value: As under 141</p> <p>See 141</p> | No Di water Error code 143 |
| 143 | <p>Adjustment of limit value fo DI water intake tvek.</p> <p>1 - 20 Min., Default = 7 Min.</p> | DI water P1 7 Minutes |
| 146 | <p>Working level with preheated DI water (M9) from boiler is not attained during fill process.</p> <p>If during intake process of preheated DI water (drain pump M32), the operating level 2 is not attained within a predetermined period of time tvew (3 min.), the control system will interrupt the current process and display error.</p> <p>Monitoring of limit value: as under 141</p> <p>Cause:</p> <p>Using control (key 3), activate K 41</p> <p>Relay K 14 is defective</p> <p>Remedy:</p> <p>replace</p> | No DI water Error code 146 |

| | | | |
|------------|---|--|--|
| | Feed pump M 14 is defective Level system leaks | repair / exchange check hoses | |
| 146 | Adjustment of limit value for preheated DI water intake t _{view} . 1 - 15 min., Default = 3 min. | | DI water boiler 3 Minutes |
| 147 | The max overflow level has been reached. If for more than seconds, the overflow level of level sensor (P1) has been exceeded, the control system will interrupt the current process and display error. (See chapter 6 Wash Chamber Level). Limit value $P_{u\ static} = 3.2\ V$ (without wash pump M 15) Limit value $P_{u\ dynamic} = P_{u\ static} - 1.3\ V$ (with running wash pump) Cause: Wash pump M 16 is deductive Pressure switch S211 stays always open Water valve doe not close properly | Remedy: Check Contact K16 or wiring connection Check plugs and wiring replace valves | Overflow level Error code 147 |
| 148 | Fill level monitoring recognizes no level change during drainage of wash chamber. If during the wash chamber drain process with drain valve Y31 or drain pump, the level has not dropped to 0.7 V (empty level + 0.2 V) within a predetermined period of time, the control system will interrupt the current process and display error. Monitoring of limit value: only during drain process until empty level is attained. Cause: The drain is plugged or drain pump is broken. The wash pump contact K 15 stays stuck closed If during the draining process (Y31 is activated) the pressure switch S211 is still closed, the control will interrupt the current process and display error. Monitoring of limit value: 5 sec. after and during drain process (5. Sec after Y31 is activated). Cause: Pressure switch S211 stays always closed Activate Y 31 (key 3) Wire break from CPU to valve Relay Y 31 does not react Drain valve defective Drain clogged (building side) Relay for wash pump K 15 remains stuck | Remedy: replace Pressure switch S211 replace control system replace clean | No drainage Error code 148 |
| 148 | Adjustment of limit value for drainage t _a . 3 - 15 min., default = 3 min. | | Drainage P1 3 Minutes |
| 149 | The minimal operating level (water heater no longer immersed) has been undercut during wash phase. Loss of water during wash phase. If during wash cycle, the minimum level of 1 V is undercut for 15 seconds, the control system will interrupt the current process and display error. Monitoring limit value: only during heating phase. Cause: | Remedy: | Level too low Error code 149 |

| | | | |
|------------|---|---|---|
| | Drain valve leaks Exhaust has excessive suction power Wtare vapor is being removed Exhaust flap remains open at all times Exhaust flap remains closed at all times Liquid is drained by means of pressure switch S31 Exhaust system is closed, expanding air cannot escape Air presses back by way of exhaust system That opens overflow level Too much air is aspirated during fill process | clean membrane or replace reduce negative pressure repair repair check exhaust system install motorized flap | |
| 150 | <p>The maximum level (water is over door threshold) has been exceeded. Undesirable water intake into wash chamber.</p> <p>If level of 3.2 V is exceeded for 3 seconds, the control system will interrupt the current process and display error.</p> <p>As long as the level is over 3.2 V, the door remains locked.</p> <p>Monitoring limit values: as long as machine is under power.</p> <p>Cause:</p> <p>Water valves are not tight</p> | <p>Remedy:</p> <p>replace water valves</p> | <p>Level too high</p> <p>Error code 150</p> |
| 151 | <p>Flow meter 1 does not recognize the required number of pulses or to high dosing.</p> <p>If in a program step phase dosing with D1 and flow meter 1 is programmed, during dosing with pump M 21, within 1 minute a minimum of 10 pulses (in the beginning it is possible that only air will come through), and after each further 6 sec. at least 50 impulses take place. The flow meter impulse is supervised during the entire dosing phase on minimum throughput: 8 Imp. per sec. of flow meter 1 must be counted. Otherwise, the control system will interrupt the current process and display error.</p> <p>The limit value for interrupt criteria X pulses in 6 seconds is calculated as follows:</p> <p>The pulse value / liter found in SW module for configuration dosage 2 used for calibrating is divided by a constant $K = 77$.</p> <p>Example:</p> <p>Pulse value / liter divided by $K = \text{limit value } X$</p> <p>Pulse value / liter = 3610</p> <p>$X = \text{limit value} = 3610 / 77 = 47$</p> <p>Error 151 is flashed if less than 47 pulses are registered within 6 seconds.</p> <p>Overdosing</p> <p>On flow meter , 6 impulses within 6 seconds must not be exceeded. Die maximum dosing is limited to 500 ml.</p> <p>Monitoring: Over full program cycle and when no dosing takes place, i.e. no output K21 – K 24 is active (Program run).</p> <p>Cause:</p> <p>Dosing module (key 2)</p> <p>Check flow meter 1</p> | <p>Remedy:</p> <p>replace flow meter 1</p> | <p>false dosing</p> <p>D1</p> |

| | | | |
|-----|--|-----------------------------|-----------------------------------|
| | Hose of dosing pump defective Dosing pump defective Dosing lines clogged IF DOSING PUMP KEY IS PRESSED, MACHINE GOES INTO ERROR MODE | replace replace clean | |
| 152 | Flow meter 2 does not recognize the required number of pulses during dosing. As per 151 | | false dosing D2 |
| 153 | Flow meter 3 does not recognize the required number of pulses during dosing. As per 151 | | false dosing D3 |
| 154 | Flow meter 4 does not recognize the required number of pulses during dosing. As per 151 | | false dosing D4 |
| 161 | <p>The temperature sensor in the wash chamber does not recognize any change in temperature during heating phase.</p> <p>If during the heating phase in the wash chamber (NTC 1), no minimal temperature change of 1°C is detected over a predetermined period of time t_h (3 min.) or during the dwell phase the temperature is no longer attained for more than t_h (3 min.), the control system will interrupt the current process and display error.</p> <p>Monitoring limit values: during heating phase with water.</p> <p>Cause:</p> <p>Activate heater contactor K 16 with control (key 3), only to be done if there is water in tank</p> <p>Heater contactor is defective</p> <p>Heater is defective</p> <p>Cable connections are burnt</p> <p>Drain valve leaks</p> <p>Exhaust system has excessive suction power</p> <p>Exhaust flap remains always open</p> <p>Remedy:</p> <p>replace</p> <p>replace</p> <p>replace connector</p> <p>clean / replace membrane</p> <p>reduce negative pressure</p> <p>repair</p> | | No heat Error code 161 |
| 161 | <p>Adjustment of limit value minimal water heating t_h.</p> <p>1 - 6 min., default = 3 min.</p> | | Delta Temp. NTC 1 3 Minutes |
| 162 | <p>The temperature sensor in DI boiler does not recognize any change in temperature during heating phase.</p> <p>If during the heating phase in the DI boiler (NTC 4), no minimal change in temperature of 1°C over a predetermined period of time t_{hve} (5 min.) is recognized, the control system will interrupt the current process and display error.</p> <p>Monitoring limit values: during heating and dwell phase.</p> | | No heat Error code 162 |
| 162 | <p>Adjustment of limit value, minimal water heating in DI boiler t_{hve}.</p> <p>1 - 6 min., default = 5 min.</p> | | Delta Temp. NTC 4 5 Minutes |
| 163 | <p>Temperature overrun</p> <p>In case water temperature overshoots target value by 8°C during water steps M 1, 2, 3, 4, 7, 9, active process is interrupted.</p> <p>Cut-off criteria:</p> | | Water too warm Error code 163 |

| | | |
|-----|--|---|
| | <p>If water temperature of NTC 1 is exceeded by 8°C for 1 minute after the fill process during a water step M 1, 2, 3, 4 or 9, the SW cuts off the current process. If no water temperature is entered (00 default), monitoring is inactive and no error is indicated.</p> <p>In function cold pre-rinse F1, max. permitted temperature is 45°C. If water temperature exceeds that value, process is also interrupted. Interrupt criteria:</p> <p>If under function F1 a temperature of 00°C is entered and water temperature of 45°C of NTC 1 is exceeded for 1 minute, the SW shuts current process down.</p> <p>Cause:</p> <p>Water media are interchanged</p> <p>Programming error of media assignments</p> <p>Heater contactor K 16 remains stuck</p> <p>Remedy:</p> <p>properly connect water</p> <p>check programs</p> <p>replace</p> | |
| 188 | <p>In configurations module 2, „Digital Interface 1“, the error display must be activated. Default = off.</p> <p>Cause:</p> <p>Interface print for Input / Output defective</p> <p>Fuse blown</p> <p>Connector has loose contact</p> <p>Coding switch in wrong position</p> <p>Remedy</p> <p>Replace interface print</p> <p>Check connector, fuse and coding switch and correct setting -> 4</p> | <p>No communication</p> <p>Error code 188</p> |
| 189 | <p>Cause:</p> <p>IPD is defective</p> <p>Door fuse blown</p> <p>Connector has loose contact</p> <p>Coding switch in wrong position</p> <p>Remedy</p> <p>Replace IPD</p> <p>Check connector, fuse and coding switch and correct setting -> 2</p> | <p>No communication</p> <p>Error code 189</p> |
| 190 | <p>Cause:</p> <p>General communication problem of SPI buss due to voltage drop, fuse door pressure print of interface print is defective, cable interruption, connector has loose contact, a code switch sits in position of another print.</p> <p>Remedy</p> <p>Fix cable break</p> <p>Replace fuse</p> <p>Check connector, fuse and coding switch and correct setting</p> | <p>No communication</p> <p>Error code 189</p> |